

- **Gold Plated BNC**
- **Pre-programmed standards**
pH: 1.68, 2.00, 4.00, 4.01, 6.87, 6.99, 9.18, 9.21, 10.01, 12.00, 12.45 (at 25°C)
Conductivity: 1413 µS/cm, 12.88 mS/cm, 111.8 mS/cm (at 25°C)
- **Stability algorithm with intuitive indicator**
- **Easy to use**
- **Ideal for schools and small laboratories**
- **Portable**
- **Accuracy at a budget**



The C5000 series are basic, yet complete and accurate portable meters. It can measure all common electrochemical parameters such as pH, ORP, Conductivity and Dissolved Oxygen. C5020 adds a larger conductivity range, Salinity, TDS and measurement storage capability.

With an optional car adaptor you can recharge batteries in the car.



Consort

SK10T



Consort

SP10T



Consort

SZ10T

Specifications depending on model

pH/MV Channels	1
EC/TDS/SAL/DO Channels	1
Temperature Channels	1
pH	0.00...14.00 pH
mV	±1000 mV
Conductivity	0...100 mS/cm 0...1000 mS/cm
Salinity	0.0...70.0
TDS	0...100 g/l
Dissolved oxygen	0.00...20.00 mg 0...200%
Temperature	0.0...+100.0°C
Warranty	36 months
Made in Belgium	

Code	Description
C1010	mV/pH/EC/DO meter
C1020	mV/pH/EC/SAL/TDS/DO meter

Code	Description	Content
C10xxP	pH meter kit, epoxy electrode	meter + pH electrode SP10T + 2x50ml pH buffers + 50ml electrolyte + carrying case
C10xxK	EC meter kit, epoxy electrode	meter + EC electrode SK10T + 50ml EC standard + carrying case
C10xxPK	pH/EC meter kit, epoxy electrodes	meter + pH electrode SP10T + EC electrode SK10T + 2x50ml pH buffers + 50ml electrolyte + 50ml EC standard + carrying case
C10xxZ	Oxygen meter kit	meter + DO electrode SZ10T + carrying case
C10xxT	Complete meter kit, epoxy electrode	meter + pH electrode SP10T + EC electrode SK10T + 2x50ml pH buffers + 50ml electrolyte + 50ml EC standard + DO electrode SZ10T + carrying case
C10xxX	Meter kit without electrodes	meter + 2x50ml pH buffers + 50ml electrolyte + 50ml EC standard + carrying case
ADAPT-CAR	Car adaptor, 12 V	

Kits are available for each meter. Replace **xx** with correct meter number. F.i. C5010P, C5020T, etc...



Artistic in science

Consort bvba
Parklaan 36
B-2300 Turnhout
Belgium

Telephone [\[+32\] \[0\]14 41 12 79](tel:+322014411279)
Email info@consort.be
Sales sales@consort.be
Support support@consort.be
Website www.consort.be